



Oncology Medical Image Database (OMI-DB): A curated platform for facilitating big data analytics in healthcare James A. Leighs¹, M. Patel¹ and M. D. Halling-Brown¹

¹Scientific Computing, Royal Surrey County Hospital, Guildford, UK

Introduction

Routine use of imaging in healthcare is continuously increasing and has been implemented in national screening programs such as Breast Cancer. Medical

Image Database

OMI-DB was designed purely to support medical research and contains anonymised unprocessed and processed images, annotations, pathology data

imaging collects a wealth of invaluable information, representing an estimated dataset of over 55 petabytes that currently remains largely unexploited. Presented here is an Oncology Medical Image Database (OMI-DB), developed as a platform for facilitating analytical healthcare research. This is an area with huge potential, particularly when placed within the context of recent cancer statistics, such as the likelihood that 1 in 2 people will be diagnosed with cancer. OMI-DB is a unique resource given the number of cases and volume of additional curated information. The data contained holds great potential for aiding machine learning studies, such as for predicting the results of invasive pathological tests.

Feature Extraction

Gabor

Laws

FeaTPy Toolkit:

Intensity Stats

A custom Python library has been written to extract first, second and higher order quantitative imaging features from images received by OMI-DB. The feature classes including Gabor, Laws and Run Length features, produce the 321 individual features summarised in the table below. The feature classes have been engineered to extract data from a variety of scales and orientations.

Run Length

and expert-determined ground truths describing regions of interest (ROIs), as shown in the example below. Currently OMI-DB holds mammographic cases but can easily be expanded to other modalities. A fully automated image retrieval system has been developed, collecting data from 4 breast screening sites in London, Cambridge, Dundee and Guildford.

> (left) an example mammogram from OMI-DB along with an expert-marked ROI

# Patients recruited	8579
# Centres recruiting	4
# Cancer cases	6358
# Cancer cases marked	3912
# Benign cases	541
# Normal Cases	1680
# Images	142,056

Texture Features:

Six feature classes extract data from images, of which two perform matrix convolution, Gabor and Law's. These higher order texture features are often used in image recognition and classification. Both construct mathematical filters, before convoluting them with the original image. The below examples show some of the different types and scales of data extracted using these features.

MedianStandard DeviationEdgesLong Run EmphasisDissimilarityVariance & RangeSkewnessSpotsGray Level Non-uniformityHomogeneitySkewness & KurtosisFrequencies = 7RipplesRun Length Non-uniformityEnergyRMSAngles = 7Additional Combinations = 6Run PercentageCorrelationEntropyImage: Standard Combinations = 6Low Gray Level Run EmphasisASMMinimumImage: Standard Combinations = 2Distances = 51st, 10th, 90th & 99th PercentileImage: Standard Combinations = 2Distances = 51414722141204	Mean	Mean	Levels	Short Run Emphasis	Contrast	Distances = 4
Variance & RangeSkewnessSpotsGray Level Non-uniformityHomogeneitySkewness & KurtosisFrequencies = 7RipplesRun Length Non-uniformityEnergyRMSAngles = 7Additional Combinations = 6Run PercentageCorrelationEntropyLow Gray Level Run EmphasisASMMinimumHigh Gray Level Run EmphasisAngles = 4MaximumDirections = 2Distances = 51 st , 10 th , 90 th & 99 th Percentile14722141204	Median	Standard Deviation	Edges	Long Run Emphasis	Dissimilarity	
Skewness & KurtosisFrequencies = 7 RipplesRipplesRun Length Non-uniformity Run PercentageEnergyRMSAngles = 7 	Variance & Range	Skewness	Spots	Gray Level Non-uniformity	Homogeneity	
RMSAngles = 7Additional Combinations = 6Run PercentageCorrelationEntropyLow Gray Level Run EmphasisASMMinimumHigh Gray Level Run EmphasisAngles = 4MaximumDirections = 2Distances = 51st, 10th, 90th & 99th Percentile14722141204	Skewness & Kurtosis	Frequencies = 7	Ripples	Run Length Non-uniformity	Energy	
EntropyLow Gray Level Run EmphasisASMMinimumHigh Gray Level Run EmphasisAngles = 4MaximumDirections = 2Distances = 51st, 10th, 90th & 99th Percentile14722141204	RMS	Angles = 7	Additional Combinations = 6	Run Percentage	Correlation	
MinimumHigh Gray Level Run EmphasisAngles = 4MaximumDirections = 2Distances = 51st, 10th, 90th & 99th Percentile14722141204	Entropy			Low Gray Level Run Emphasis	ASM	
Maximum Directions = 2 Distances = 5 1 st , 10 th , 90 th & 99 th Percentile 147 22 14 120 4	Minimum			High Gray Level Run Emphasis	Angles = 4	
1st, 10th, 90th & 99th Percentile 14 147 22 14 120 4	Maximum			Directions = 2	Distances = 5	
14 147 22 14 120 4	1 st , 10 th , 90 th & 99 th Percentile					
	14	147	22	14	120	4

Pipeline:

- The custom feature extraction toolkit has been implented within a data extraction pipeline for automated and ondemand feature extraction. When new data is received by
- OMI-DB, the feature pipeline is



GLCM

LBP



Conclusions

The main focus has been on developing automated image collection and feature extraction, resulting in a large, curated database of images and associated feature data. Preliminary efforts to analyse potential for predicting results of pathological tests have showed promise. Law's and Gabor features have shown to be the most significant biomarkers. Two diagnostic/prognostic indicators have been investigated. Disease grade, a measure of the rate at which the cancer is likely to spread, has been predicted with accuracies up to 65%. Molecular subtype, a diagnostic classifier used to target personalised treatment plans, has been predicted up to 82% accuracy. Further work is ongoing, aiming to both improve classification rates and predict further pathological results. Success in this area has potential benefits in prioritsing and streamlining patient care by reducing wait times for results. Impacts of this could be significant, given early diagnosis is seen as one of the most important factors for survival.

run and the full set of image features are extracted from each image and corresponding

lesion ROI.

Feat-DB:

(above) a schematic of the automated feature extraction pipeline (1) image details received (2) image data is transferred from OMI-DB (3) ROI details are queried from a local clinical database (4) FeaTPy feature extraction toolkit is called (5) feature data is returned (6) features are written to Feat-DB or local CSV file

FeaTPY currently extracts 321 features from each image and ROI, including an add-on to adjust ROI size and accept CAD segmentation co-ordinates. There are >100,000,000 features contained within Feat-DB., including quantitative imaging features, additional pathological data and other annotations.

Acknowledgements

OMI-DB development was undertaken as part of the OPTIMAM project and supported by CR-UK

Corresponding Author Email: mishalpatel@nhs.net